


# Mammalian expression vector pHcRed-Tandem-N1



For vector sequence, please visit our Web site at www.evrogen.com/support/vector-info.shtml

Product	Cat.#	Size						
pHcRed-Tandem-N1	FP204	20 µg						

Please contact your local distributor for exact prices and delivery information.

Vector type mammalian expression vector

Reporter t-HcRed
Reporter codon usage mammalian
Promoter for t-HcRed P<sub>CMV IE</sub>
Host cells mammalian

Selection prokaryotic — kanamycin

eukaryotic — neomycin (G418)

Replication prokaryotic — pUC ori eukaryotic — SV40 ori

## Multiple cloning site (MCS)

Nhe I	Bgl II*	Sac I	Hind III	EcoR I		Sal I		Kpn I	_	Apa I		BamH	<u> </u>	Age I			T-HCR	<b>=</b>	
GCTA.GCG.CTA.CCG.GAC.	TCA.GAT.CTC.	GAG.CTC.	AAG.CTT.	CGA.ATT.	CTG. CAG	.TCG.	ACG.	GTA.C	CG.CGG	.GCC.	CGG.	GAT.	CCA.C	CCG.GTC	. GCC	ACC.	ATG. GTG	. AGC	
Eco47 III	Xho	1*			Pst I				Sac II	Sma I	* /Xm	— a I*							

<sup>\* -</sup> not unique sites.

### Use

- Generation of fusions to the t-HcRed N-terminus
- Expression of t-HcRed or its fusions in mammalian cells

# **Vector description**

pHcRed-Tandem-N1 is an eukaryotic (mammalian) expression vector encoding far-red fluorescent protein t-HcRed. The vector is designed to generate t-HcRed-tagged fusions or to express t-HcRed in mammalian cells.

t-HcRed codon usage is optimized for high expression in mammalian cells (humanized, Haas et al., 1996). To increase t-HcRed translation, Kozak consensus translation initiation site is generated upstream of t-HcRed sequence (Kozak, 1987). Multiple cloning site (MCS) is located between  $P_{\text{CMV IE}}$  and t-HcRed coding sequence.

The vector backbone comprises immediate early promoter of cytomegalovirus (P<sub>CMV IE</sub>) for protein expression, SV40 origin for replication in mammalian cells expressing SV40 T-antigen, pUC origin of replication for propagation in *E. coli* and f1 origin for single-stranded DNA production. SV40 polyadenylation signals (SV40 poly A) direct proper processing of the 3' end of the reporter mRNA.

SV40 early promoter provides neomycin resistance gene expression to select stably transfected eukaryotic cells using G418. Bacterial promoter (P) provides kanamycin resistance gene expression in *E. coli*. Kan<sup>r</sup>/Neo<sup>r</sup> gene is linked with herpes simplex virus (HSV) thymidine kinase (TK) polyadenylation signals.

# Generation of fusions

A localization signal or a gene of interest should be cloned into MCS of the vector. It will be expressed as a fusion to the t-HcRed N-terminus when inserted in the same reading frame as t-HcRed and no in-frame stop codons are present. The inserted sequence should contain an initiating ATG codon. t-HcRed-tagged fusions retain fluorescent properties of the native protein allowing fusion localization *in vivo*.

**Note:** The plasmid DNA was isolated from dam<sup>+</sup>-methylated *E.coli*. Therefore some restriction sites are bloked by methylation. If you wish to digest the vector using such sites you will need to transform the vector into a dam<sup>-</sup> host and make fresh DNA.

#### **Expression in mammalian cells**

pHcRed-Tandem-N1 can be transfected into mammalian cells by any known transfection method. If required, stable transformants can be selected using G418 (Gorman, 1985).

Unmodified pHcRed-Tandem-N1 will express t-HcRed, when transfected into eukaryotic (mammalian) cells.

## Propagation in E. coli

Suitable host strains for propagation in *E. coli* include DH5alpha, HB101, XL1-Blue, and other general purpose strains. Plasmid incompatibility group is pMB1/ColE1. The vector confers resistance to kanamycin (30 µg/ml) to *E. coli* hosts. Copy number in *E. coli* is about 500.

#### **Location of features**

P<sub>CMV IE</sub>: 1-589

Enhancer region: 59-465 TATA box: 554-560

Transcription start point: 583

MCS: 591-671 t-HcRed gene

Kozak consensus translation initiation site: 672-682

Start codon (ATG): 679-681; Stop codon: 2086-2088

Polyadenylation signals: 2218-2223 & 2247-2252

mRNA 3' ends: 2256 & 2268

f1 single-strand DNA origin: 2315-2770 (packages the noncoding strand of t-HcRed)

Bacterial promoter expression of Kan<sup>r</sup> gene:
-35 region: 2832-2837; -10 region: 2855-2860

Transcription start point: 2867

SV40 origin of replication: 3111-3246 SV40 early promoter

Enhancer (72-bp tandem repeats): 2944-3015 & 3016-3087 21-bp repeats: 3091-3111, 3112-3132 & 3134-3154

Early promoter element: 3167-3173

Major transcription start points: 3163, 3201, 3207 & 3212

Kanamycin/neomycin resistance gene

Neomycin phosphotransferase coding sequences: Start codon (ATG): 3295-3297; stop codon: 4087-4089

G->A mutation to remove PstI site: 3477

C->A (Arg to Ser) mutation to remove BssHII site: 3823

Herpes simplex virus (HSV) thymidine kinase (TK) polyadenylation signal

Polyadenylation signals: 4325-4330 & 4338-4343 pUC plasmid replication origin: 4674-5317

# References

Gorman C. (1985) In DNA cloning: A Practical Approach, Vol. II. Ed. D. M. Glover. (IRL Press, Oxford, U.K.), pp. 143-190.

Haas J. et al. (1996) Curr. Biol. 6: 315-324.

Kozak M. (1987) Nucleic Acids Res. 15:8125-8148.

### Notice to Purchaser:

Evrogen Fluorescent Protein Products (the Products) are available to Purchasers for non-commercial non-for-profit research use. With purchase of the Products, Purchaser is granted a worldwide, non-exclusive, royalty-free, limited license to use the Products for non-commercial life science research only. Such license specifically excludes the right to sell or otherwise transfer the Products, its components or derivatives to third parties and any uses or activities (or the results therefrom) that themselves generate revenue for the Purchaser. For commercial use of the Products please contact Evrogen at *license@evrogen.com* for license information.

The CMV promoter is covered under U.S. Patents 5,168,062 and 5,385,839 and its use is permitted for research purposes only. Any other use of the CMV promoter requires a license from the University of Iowa Research Foundation, 214 Technology Innovation Center, Iowa City, IA 52242.

# MATERIAL SAFETY DATA SHEET INFORMATION

To the best of our knowledge, these products do not require a Material Safety Data Sheet. However, all the properties of these products (and, if applicable, each of their components) have not been thoroughly investigated. Therefore, we recommend that you use gloves and eye protection, and wear a laboratory coat when working with these products.